
QMultiTractor
Type Karlsruhe

A New Way of Alveolar Distraction

High stabilisation by minimally invasive technique

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This product was designed with friendly support of **Oral and Maxillofacial Clinic in Karlsruhe**

Plastic-, reconstructive, and esthetic Surgery - Dental Ambulancy
Director Priv.-Doz. Dr. Dr.(H) A. Dunsche

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
The most important characteristics of Q-MultiTractor, Type Karlsruhe

- Solid, modular mini plate design, comprising distraction spindle (A), transport plate (B) and base plate (C)
- many combinations possible fitting to needed indication
- 20 mm distractionheight possible with the 10 mm spindle
- modularity enables considerable easy and safe handling
- minimally invasive surgery, especially when detached
- outstanding stability against lingual movement while using the innovative pin-baseplate
- apparatus made of titanium- thereby providing minimum potential for allergisation
- neatly defined instrumentarium

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
Indications

Distraction osteogenesis resembles a modern method for the absolute elevation of the alveolar process, thus also presenting an alternative solution to the conventional augmentation treatment of the alveolar process. Contrasting to the augmentation via bone substituting or bone building materials or osseous grafts, this technique makes use of the on-site bone following segmental osteotomy, the device then expanding the bony tissue and forming a callus which is stimulated for growth. The distraction procedure does not only dilate the osteomised bone around the osteotomic fissure, but serves soft tissue stretching as well. This desirable phenomenon may be utilised not only for dental implantology, but also for preprosthetic and reconstructive surgery.

Contraindications

- Height of the alveolar process below 8mm
- Proximity to important anatomic structures (with a spacing of less than 2 mm, e.g. to the inferior dental nerve)
- Severe osteoporosis
- General contraindications due to severe system diseases

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Patient key:


- small size helps within difficulty accessible ranges
- drilling hole provides aspiration protection
- lasermarked arrow indicates the right direction of rotation
- one full turn makes 0,35 mm distracted height

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

- High-precision conical thread**
- avoids rotation movements
 - fixes distractor spindle

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- Pin to stabilize distractionvector**
- neutralises lingualisation
 - allows for minimally invasive surgery

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
20 mm Distractionheight with 10 mm Spindle

The distraction phase is followed by the retention phase. To make this process as comfortable as possible for the patient, the distraction spindle is countersunk. To do so, a stab incision is used to expose and loosen the fixation screw on the transport plate. It is then possible to return the distraction spindle to its original position by turning against the distraction direction. Next the transport plate is tightened via the fixation screw.

Only modular systems like the **Q-MultiTractor, Type Karlsruhe**, allow for this form of countersinking.

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abb. / fig.7

abb. / fig.8

The picture sequence demonstrates the use of Q-MultiTractor with single plate

This model operation starts with a groove near the symphysis of corticalis, mesial the exit of nervus mentalis (Fig. 1).

The hole for the pinplate is carried out from vestibulum to lingual, using the drilling gauge to align the distraction direction (Fig. 2).

Then the pin plate can be inserted easily and fixed with two screws. (Fig. 3).

The distraction spindle is screwed inside the pinplate, the transport plate adjusted onto the spindle, adapted and fixed with 3-4 screws. As a general rule the first pierced hole of the transport plate is left blank to facilitate the adjustment of the plate. (Fig. 4).

Contrary to the compact distractor only the distractor spindle is unscrewed now temporarily. Horizontal bone cutting takes place after removal of the spindle. Therefore the transport plate is used as guidance. For vertical cutting it is favorable to carry out first only one of the two cuttings. (Fig. 5,6).

After refixing transport plate and spindle the second vertical bonecut is set. This procedure assures that the transport segment is always partially fixed. Both vertical cuttings should form an open angle. (Fig. 7).


With final test distraction the model operation should be reviewed for mobility (Fig. 8).

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
The picture sequence demonstrates the use of Q-MultiTractor with double plate

This model operation starts with a groove near the symphysis of corticalis, mesial the exit of nervus mentalis (Fig. 1).

The hole for the pinplate is carried out from vestibulum to lingual, using the drilling gauge to align the distraction direction (Fig. 2).

The distraction spindle is screwed inside the pinplate, the transport plate adjusted onto the spindle, adapted and fixed with 3-4 screws. As a general rule the first pierced hole of the transport plate is left blank to facilitate the adjustment of the plate. (Fig. 3).

Contrary to the compact distractor only the distractor spindle is unscrewed now temporarily. Horizontal bone cutting takes place after removal of the spindle. Therefore the transport plate is used as guidance. For vertical cutting it is favorable to carry out first only one of the two cuttings. (Fig. 4,5,6,7).

After refixing transport plate and spindle the second vertical bonecut is set. This procedure assures that the transport segment is always partially fixed. Both vertical cuttings should form an open angle. (Fig. 7).

With final test distraction the model operation should be reviewed for mobility (Fig. 8).

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

1. Clinical case:

Starting situation: with loss of teeth and alveolar bone (fig. 1). The operation begins with incision of the mucous membrane in the vestibulum. The bone is exposed afterwards.

The placement of the basis plate is prepared with the rotating drill. The bone drilling for the admission of the pin basis plate takes place with the help of the drilling gauge. (fig. 2)

The basis plate with the distraction spindle is inserted into the bone and fixed. (fig. 3)

After adaptation of the transportation plate, it is fixed with 3-4 screws per segment. Then the distraction spindle is removed with the help of the holding key. (fig. 4)

In the next step the horizontal bone cut takes place with a saw (fig. 5). The transport plate serves as guide. Afterwards, the distraction spindle is reinserted and fixed in transport and basis plate.

After the final vertical cuts are executed in an open angle the test distraction takes place to assure the all-round transection of the bone. (fig. 6)

The operation ends with a double layer suturing of mucous membrane. (fig. 7)

Figure 8 shows the situation at the end of the distraction phase (18 days post OP).

During the distraction phase the distraction spindle lifts itself (fig. 8, 9) and is sunk in the retention phase again. (fig. 10)

In a second operation the Q-MultiTractor is completely removed.

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



2. Clinical case:

Starting situation: strongly atrophic lower jaw, loss of teeth and alveolar bone. (Fig. 1)

Very small bone height prevents classical proceeding of the distraction. With the existing width of the bone only the possibility of a visor osteotomy using the Q-MultiTractor is given.

The operation is prepared in detail with radiographs (fig. 2) and a 3D-reconstruction. (fig. 3, 4)

After the bone is exposed, a wide basis plate is fixed with several screws.

Subsequently, the transport plate is placed onto the distractions pindle and bent approx. 90° to lingual and fixed to the bone. The distraction spindle is now removed.

With the perpendicular, arc-shaped cut, following the form of the lower jaw, the visor osteotomy is accomplished. (fig. 5, 6)

After reinserting the distraction spindle the test distraction takes place to secure a perfect mobility of the transport segment. (fig. 7)

After the double layer suturing of mucous membrane the operation is completed. (fig. 8)

Subsequently, x-ray control takes place, the correct positioning of the Q-MultiTractor and the distraction vector can be examined. (fig. 9)

Further x-ray controls during the distraction phase monitor the progress. (fig. 10)

At the end of the distraction sufficient new bone is created, the Q-MultiTractor can be removed. (fig. 11)

After removal of all distractor elements the new bone is visible. (fig. 12)

Afterwards, an implantation took place with four QK-implants in this case. (fig. 13, 14)

q-Implant® System

- Q-Implant-System fits all indications
- ideal for immediate loading and temporisation
- high primary stability, because of the special thread design
- two-phase implants are compatible to many other systems

q-bone-grafting-set

- suitable for onlay-plastic and mesh supported augmentation
- bone screws in diameter 1,0 und 1,3 mm
- colour coding
- including titanium bowl for blending augmentation material

q-mesh

- 3-dimensional formed titanium mesh
- applicable to totally atrophied maxilla
- reduced operation time
- easy to use
- individual modulation

bone-pin-system

- titanium pins in length 3 mm and 5 mm
- for attachment of titanium mesh, -foils and membranes
- appropriate for 3-dimensional bone reconstruction
- titanium mesh 0,1 und 0,2 mm
- titanium foil 20 µ und 40 µ

osteosynthesis system

- appropriate for maxillofacial-surgery
- titanium bone screws and bone plates
- screw diameter 1,0/1,2/1,5/1,8/2,0 and 2,3 mm
- screw top optional with inner cross head or inner square
- large variety of plates with thickness 0,6 mm and 1,0 mm

scalpels

- scalpel blades
- disposable scalpels

multi-f

- external Ortho-Fixator made of titanium
- modular component system
- easy to use
- distraction, compression, dynamisation

HIGH QUALITY TITANIUM PRODUCTS SINCE 1993

MEDICINE

INDUSTRY

CONSUMER

TRINON
TITANIUM

TRINON Titanium GmbH · Augartenstraße 1 · D-76137 Karlsruhe
Tel.: +49 721 93 27 00 · Fax: +49 721 24 991 · www.trinon.com