

1




A Guide Drill (Diameter 2mm) is used to drill pilot holes into the plastic model.

2



Guide Pins (Diameter 2 mm) are inserted into the guide holes to check drill depth and to locate the Master Tube and Drill Guides.

3



The Master Tube and Drill Guide (Diameter 2 mm) are inserted into the mold.

4



To fix the Master Tube with a surgical stent, resin is applied onto the model.

5



Finally the surgical stent is removed from the plastic model and the Drill Guides are replaced as required by the implant system. If required, Manual Drill Guides can be used with the RGIT.

APPLICABLE FOR ALL IMPLANT SYSTEMS

The Reverse Guide Implant Technique utilizes surgical guides without the use of any software.

It offers surgeons an easy preparation, high precision and increased safety of implant placement.

The technique is practice based and offers a fast and low-cost solution for doctors.

Only one disposable Master Tube is required per implant.

The technique is easy to correct making it ideal for beginners and an invaluable aid for professionals

For more information please contact:

TRINON
TITANIUM

Trinon Titanium GmbH
Augartenstr. 1
D-76137 Karlsruhe / Germany

Tel.: +49 721 932700
Fax: +49 721 24991

trinon@trinon.com
www.trinon.com

Reverse Guide Implant Technique



TRINON
COLLEGIUM
PRACTICUM

Product

Art. No.

Product Name



Q6001

RGIT Drill Diameter 2mm



Q6002

RGIT Pin Diameter 2mm



Q6003

RGIT Master Tube Diameter 4mm



Q6020

RGIT Drill Guide Diameter 2mm

Q6022

RGIT Drill Guide Diameter 2.2mm

Q6025

RGIT Drill Guide Diameter 2.5mm

Q6028

RGIT Drill Guide Diameter 2.8mm

Q6030

RGIT Drill Guide Diameter 3.0mm

Q6033

RGIT Drill Guide Diameter 3.3mm

Q6035

RGIT Drill Guide Diameter 3.5mm



Q6120

RGIT Drill Guide Manual Diameter 2mm

Q6122

RGIT Drill Guide Manual Diameter 2.2mm

Q6125

RGIT Drill Guide Manual Diameter 2.5mm

Q6128

RGIT Drill Guide Manual Diameter 2.8mm

Q6130

RGIT Drill Guide Manual Diameter 3.0mm

Q6133

RGIT Drill Guide Manual Diameter 3.3mm

Q6135

RGIT Drill Guide Manual Diameter 3.5mm