

Introducing the new Self Drilling Twisted Trocar Screw



Twisted Trocar detail

No more drilling

With its new geometry designed by Trinon Titanium the Twisted Trocar Self Drilling Screw is ideally suited as replacement for applications where drilling is required.

The TT screw is easily placed in any bone, even cortical bone.

How it's done

This new screw type is based on the trocar principle with an added geometrical feature to enhance insertion and reduce the strength needed for fixation. With a specialized twist of the trocar tip the TT screw offers surgeons a new tool for fixation into any bone type with ease.

The TT bone screws are available in any size, any head configuration and any thread pitch from Trinon Titanium/Germany.



TRINON
TITANIUM